



MOMO

ML-BASED
ALPHA TRACK ANALYZER
BY: PANNATHAD


HANDBOOK

Software Details:	MODEL: Ultralytics YOLOv8	Measured properties Total counts (unit: count) Area (unit: pixel and μm^2) Equivalent_diameter (unit: pixel and μm^2) Mean_intensity* (unit: none) solidity* (unit: none)
	Detecting speed: 275.44 - 350 ms per tilted image. 5.5088 - 7 second per image. (For image with resolution 2992x2000)	
	Output Type: .CSV .TXT .JPEG	
	Version: Alpha v2.0.1 (2023 OCT 24)	
Contact:	peak.panthanuwong@gmail.com	Default Setting iou = 0.1 confident level = 0.21 scale = 0.6 micron/pixel tiled image size = 640x640

- This software is in the development process; any error or stability issue might occur. If there are problems or errors, please send the details to: peak.panthanuwong@gmail.com.

USER-INTERFACE



- Indicate whether they have been resolved or are still pending action



HOW TO USE

STEPS		Status	Notes
1: Select your scale ratio and additional settings (optional).		OPTIONAL	How to find a scale ratio guide is located below in the next section.
2: Choose where you want to save your result.		MANDATORY	The result will be saved in the format 'date-time.csv'.
3: Locate the image's location and start the detection process.		START	The selected folder should contain only .JPEG file
4: Wait until the program finishes analyzing.		COMPLETED	The result will be saved at the location you selected in the first step.

FINDING SCALE RATIO

STEPS	Task Owner
	<p>1.Take a scale image through the microscope.</p>
	<p>2.Run a program and use this image as the one you want to detect.</p>
 <p>sized: 640x640 pixels</p>	<p>3. After the program has run completely, go to the tiled images folder and pick the image in full scale.</p>
	<p>4.Measure the length (micron) of the image and divide by 640.</p>

RESULT FILE

- **Date-Time.csv**

Show the properties of each detected track.

	A	B	C	D	E	F	G	H
1	Name	Count	area	equivalent_diameter	mean_intensity	solidity	area_sq_microns	equivalent_diameter_microns
2	DSC_0051_0_0.JPG	1	895	33.75721245	237.0502793	0.983516484	322.2	20.25432747
3	DSC_0051_0_1000.JPG	3	2105	51.77035099	242.7648456	0.986872949	757.8	31.06221059
4	DSC_0051_0_1000.JPG	3	1805	47.93951792	241.5789474	0.987958402	649.8	28.76371075
5	DSC_0051_0_1000.JPG	3	1062	36.77200561	237.2316384	0.985157699	382.32	22.06320337
6	DSC_0051_0_1500.JPG	2	1305	40.76245338	239.3678161	0.988636364	469.8	24.45747203
7	DSC_0051_0_1500.JPG	2	198	15.87770229	216.3636364	0.975369458	71.28	9.526621371
8	DSC_0051_0_2000.JPG	2	685	29.53250901	234.5255474	0.981375358	246.6	17.7195054
9	DSC_0051_0_2000.JPG	2	966	35.07063444	239.1614907	0.991786448	347.76	21.04238067
10	DSC_0051_0_500.JPG	3	796	31.83549399	235.138191	0.987593052	286.56	19.10129639
11	DSC_0051_0_500.JPG	3	2080	51.46200786	242.7403846	0.99047619	748.8	30.87720472
12	DSC_0051_0_500.JPG	3	2063	51.25127492	242.886573	0.991826923	742.68	30.75076495
13	DSC_0051_1000_1000.JPG	2	581	27.19838553	232.1772806	0.976470588	209.16	16.31903132
14	DSC_0051_1000_1000.JPG	2	1539	44.26641683	240.4191033	0.983386581	554.04	26.5598501
15	DSC_0051_1000_1500.JPG	3	1716	46.7426899	241.7744755	0.984509466	617.76	28.04561394
16	DSC_0051_1000_1500.JPG	3	1738	47.04136827	243.1156502	0.988623436	625.68	28.22482096
17	DSC_0051_1000_1500.JPG	3	322	20.24804024	223.3229814	0.981707317	115.92	12.14882414
18	DSC_0051_1000_2000.JPG	6	2080	51.46200786	242.6177885	0.989533777	748.8	30.87720472
19	DSC_0051_1000_2000.JPG	6	1299	40.66863864	239.8845266	0.9878327	467.64	24.40118318
20	DSC_0051_1000_2000.JPG	6	418	23.06976657	226.3277512	0.976635514	150.48	13.84185994
21	DSC_0051_1000_2000.JPG	6	1243	39.78236738	237.767498	0.987291501	447.48	23.86942043
22	DSC_0051_1000_2000.JPG	6	1103	37.47510131	239.5104261	0.977836879	397.08	22.48506078

- **Date-Time_result.csv**

Show total counts and average values for each image

	A	B	C	D	E	F	G	H
1	Name	Total_count	Avr Area	Avr diameter	Avr intensity	Avr solidity	Avr area in micron	Avr diameter in micron
2	DSC_0055.JPG	44	1372.568182	40.71026682	239.3105922	0.983846625	494.1245455	24.42616009
3	Name	Total_count	Avr Area	Avr diameter	Avr intensity	Avr solidity	Avr area in micron	Avr diameter in micron
4	DSC_0056.JPG	50	1222.84	38.33991619	238.2002796	0.985062078	440.2224	23.00394971
5	Name	Total_count	Avr Area	Avr diameter	Avr intensity	Avr solidity	Avr area in micron	Avr diameter in micron
6	DSC_0057.JPG	42	1210.785714	38.34002398	238.2985947	0.985296445	435.8828571	23.00401439
7	Name	Total_count	Avr Area	Avr diameter	Avr intensity	Avr solidity	Avr area in micron	Avr diameter in micron
8	DSC_0058.JPG	52	1494.403846	42.51343038	239.4202603	0.984639986	537.9853846	25.50805823
9	Name	Total_count	Avr Area	Avr diameter	Avr intensity	Avr solidity	Avr area in micron	Avr diameter in micron
10	DSC_0059.JPG	51	1332.392157	40.2971705	239.9719434	0.983343718	479.6611765	24.1783023
11	Name	Total_count	Avr Area	Avr diameter	Avr intensity	Avr solidity	Avr area in micron	Avr diameter in micron
12	DSC_0060.JPG	58	1372.258621	41.03589551	240.1124479	0.985884918	494.0131034	24.62153731

- **Predict images (.JPEG)**

Located in “MOMOAnalyzer\runs\segment\predict“

Show images with detected tracks highlighted in colors.



• .txt

Show track type and coordinate of segmented area in tensor format.

Located in “MOMOAnalyzer\runs\segment\predict\labels”